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Ⓕ Mode d'emploi

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ⒼⒷ Instructions for use 

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**boso**  
medicus

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---

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## Equipment and documentation supplied

---

- 1 Sphygmomanometer  
**boso-medicus**
- 1 Cuff CA 01
- 1 Protective case
- 1 Blood pressure log book
- 1 Guarantee card
- 4 Batteries round cell LR6
- 1 Instructions for use

## Introductory information

---

Your **boso-medicus** is a sphygmomanometer computer that has been designed and constructed on the basis of the latest technological developments. It operates on the **oscillometric principle**. The pressure variations (oscillations), resulting from the pulse waves and transmitted by the cuff, are stored and evaluated by a microprocessor. The major advantage of this method is that no microphone is required; where as normally the reliability of a blood pressure measurement is strongly dependent on the accuracy with which the microphone is positioned when the cuff is put in to place.

This sphygmomanometer complies with the European directives underlying legislation on medical products (CE symbol), as well as the European Standard EN 1060 Part 1: 'Non-invasive Sphygmomanometers - General Requirements' and Part 3: 'Supplementary Requirements Governing Electronic Blood Pressure Measurement Systems'.

Regular checks on the accuracy of the measurement system should be carried out in accordance with the valid directives of the individual countries, governing this issue.

The technical testing procedures are detailed on page 60.

## Blood pressure readings

**To determine correct blood pressure, the following two values have to be ascertained:**

- The **Systolic or upper blood pressure**. This is created when the heart muscle contracts and blood is being forced into the blood vessels.
- The **Diastolic or lower blood pressure**. This is created when the heart muscle extends and once again fills with blood.

Measurement values are expressed in **mmHg (mm mercury column)**.

The World Health Organization (WHO) has set the following guideline for the assessing of blood pressure values:

	Systolic	Diastolic
Too high	over 140 mmHg	over 90 mmHg
Normal-borderline	130 to 139 mmHg	85 to 89 mmHg
Normal	120 to 129 mmHg	80 to 84 mmHg
Optimum	up to 119 mmHg	up to 79 mmHg


Whether medical treatment is necessary depends not only on the blood pressure but also on the risk profile of the individual patient. Please contact your GP if one of the values (SYS/DIA) constantly exceeds the limit (> 140/> 90).

## Preliminary remarks on measuring your own blood pressure

---

- The measurement should always be performed when at rest. It is recommended that the blood pressure is measured twice a day - in the morning after getting up, and in the evening when you have relaxed after a day's work.
- Blood pressure readings should be performed while sitting. Sit down in a comfortable chair and relax. The chair should have arm rests or be positioned at a table so that your arm can be supported in a slightly bent position.
- The measurement is carried out on the bare upper arm. It is important to ensure that the arm is not constricted by tight fitting clothing that has been pushed back or rolled up. If necessary remove the article of clothing before starting to measure your blood pressure.
- During the measurement, the arm with the cuff must not be moved.
- An interval of at least two minutes should be observed, before repeating a measurement.
- Variations in blood pressure are perfectly normal. Even a repeat measurement can produce considerable differences. One-off or irregular measurements are therefore of little value when trying to assess the actual blood pressure. Reliable assessment of blood pressure is only possible if you monitor your blood pressure regularly under similar conditions and record the readings in your blood pressure log book.
- Measuring your own blood pressure is not, however, a therapy. Never alter of your own accord the dosage of any medication that your doctor has prescribed for you.
- Cardiac dysrhythmia (irregular heartbeat) can reduce the accuracy of, or even give wrong readings. Measurements may also be impaired in the case of those patients with a pacemaker who have a weak pulse. The sphygmomanometer does not have any influence on the pacemaker.
- Faulty measurements can also occur if the unit is operated in the vicinity of strong electromagnetic fields (e.g. radiological applicators, mobile phones, etc.)
- Medical Electrical Equipment needs special precautions regarding EMC. You will get appropriate information from the supplier on demand.

## Setting-up the unit

### Inserting the batteries

The battery compartment is located on the underside of the unit. Insert the batteries in the manner shown in the picture above, ensuring correct polarity.

Warning: If the batteries are inserted incorrectly, the unit will not function. This can cause heat to be generated in the battery case and may cause the batteries to leak.

- Only use high-quality leak-proof batteries of the specified type (see 'Technical data'). Never mix old and new batteries, or batteries from different manufacturers.
- Remove the batteries if the unit is not going to be used for some time.
- The batteries require replacing when the symbol  is illuminated in the display.
- **Protect the environment!**  
Empty batteries and storage batteries should not be thrown away with household waste. They can be handed in at a battery collection point or disposed of as special waste. Please contact your local authority for details.

## Preparing to measure blood pressure

### How to wear the cuff

- Open cuff into a ring shape and pull onto your upper arm until the lower cuff edge lies 2-3 cm (1 inch) above the elbow. The cuff is positioned correctly if the special marking is located directly on the artery. The metal ring must never come to lie on the artery as this may cause a false reading.

Most persons have a higher blood pressure on their left arm. Therefore the measurement should be taken on the left arm. If the blood pressure is higher on the right arm, the measurement should be taken on the right arm.

- Pull the end of the cuff, which has been fed through the metal bracket, and wrap it round your upper arm. By pressing down lightly the cuff will be secured in position by the Velcro fastener.

The cuff should not be too tight; it should be possible to place two fingers between your arm and the cuff.

- Note on cuff size:


The unit must only be used in combination with the following cuffs:

Model	Rubber bag size	Circumference	Order number
CA01	13 x 25 cm	22 – 32 cm	143-4-750*
CA02	14 x 30 cm	32 – 48 cm	143-4-757

\*supplied with the unit

### Important information:

A blood pressure measurement should not take longer than 2 minutes. If a problem occurs you can remove the cuff at any time by simply releasing the Velcro fastening.


## A step-to-step guide to measuring your blood pressure

- Attach the cuff to the monitor by pressing the plug (click) into the socket on the left side of the unit.


- Set the pressure selector switch '**DRUCKVORWAHL**' to the desired cuff pressure which should be at least 30 to 40 mmHg higher than the systolic (upper) blood pressure reading. If you are not sure of this value, set the selector switch to 150 mmHg for the first measurement.

If this pressure is insufficient, the unit will switch over to the next higher pressure setting automatically.

- With the cuff in position, place your arm on the table and relax. Bend your arm slightly so that the cuff is at the same height as your heart.
- Press the START/STOP button. All display symbols light up briefly to test correct functioning.

- The unit now inflates the cuff to the preselected cuff pressure. The rise in pressure in the cuff is indicated on the digital display.
- An acoustic signal sounds once the preselected pressure has been reached. The pump switches off and the air in the cuff escapes automatically.

- During the measurement phase, the falling pressure and the symbol  are indicated on the display.
- If not already the case, you should now keep your arm completely still and not speak.

## Reading the display

- As soon as the unit registers the first beats, the  symbol blinks in time with the pulse rhythm and an acoustic signal can be heard.
- The end of the measurement is indicated by a long beeping sound. At the same time, the built-in solenoid valve opens automatically to rapidly release the remaining pressure in the cuff.
- The measurement readings shown in the display.

Systole mmHg

Diastole mmHg

Pulse/min.

- The unit switches off automatically after approximately one minute. To repeat a measurement, press the START/STOP button again.
- If you do not want to take a further measurement, remove the cuff from the arm.
- Enter the readings into the blood pressure log book. The readings will be displayed for one minute. Then the unit switches off automatically.

### Aborting a measurement:

A blood pressure measurement should not take longer than two minutes. A measurement can be aborted at any time by simply pressing the START/STOP button which automatically deflates the cuff. (To perform another measurement press the START/STOP button again.)

## Mains operation

---

- A mains connector socket is located on the rear side of the unit. If you wish to operate the unit from the mains, use only the **boso-power supply pack (order number: 410-7-150)**. This power supply pack produces a stabilized output of the correct polarity specially designed for use with the **boso-medicus** unit. Other commercially available power supply units may damage the electronic components and lead to the guarantee becoming void.


Please remove the batteries before operating the unit from the power supply pack.

Occasionally faults in the mains network may cause too high a pulse rate to be indicated. If the beeping sound becomes irregular despite the fact that the pulse rate is steady, the measurement should be repeated (if problems persist, switch to battery mode).

## Error indicator

- If a problem should arise during a measurement which prevents a usable measurement from being made, an error number is displayed in place of the usual measurement result.

- Error number key:

Err 1: Inflation pressure inadequate after three attempts to inflate.

→ Select the next higher pressure setting.

Err 2: No pulses were detected that could be evaluated.

→ Check the positioning of the cuff.

Err 3: Leak in the air system.

→ Check air connector on cuff. Cuff possibly defective.

## Cleaning and care of the unit

- Please use only a soft, dry cloth to clean the unit.
- Small stains on the cuff can be carefully removed with a stain remover.

## Guarantee information, Care of the Environment

For details on the conditions of guarantee, please refer to the guarantee card supplied.

A guarantee claim can only be accepted if the card, which has been filled-in and bears the stamp of the place of purchase, is sent in together with the unit.

For both guarantee and repair work, please send the instrument carefully packed and postage paid either to your authorized dealer or directly to:

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Used batteries and defective electronic devices do not belong into the household waste.

They can be deposited at special centres. Please contact your local authorities for information.

## Technical data

---

**Measurement principle:** oscillometric

**Measurement range:** 40 to 280 mmHg  
40 to 200 pulses/min

**Cuff pressure:** 0 to 320 mmHg

**Display:** LCD

**Operating conditions:**

ambient temperature: 10°C - 40°C

relative humidity: 10% - 85%

**Transport and storage conditions:**

temperature: -5°C - +50°C

relative humidity: max. 85%

**Power supply:**

DC 6V (4 x 1.5 V round cell IEC LR 6, alkaline manganese)

Alternatively as an extra: power supply pack DC 6V,

order no.: 410-7-150

(polarity: outer MINUS, inner PLUS) 

**Battery check:** indicator symbol in display

**Weight:** 0.6 kg

**Classification:** Protection class II (回)

Type BF (人)

**Clinical trials:** DIN 58130 testing procedure N2  
'Comparative Blood Pressure  
Measurements on the Left and Right Arm'

**Results:**

● systematic measurements variance:

systolic blood pressure: -0.98 mmHg

diastolic blood pressure: -3.2 mmHg

● empirical standards variance:

systolic blood pressure:  $\pm 5.96$  mmHg

diastolic blood pressure:  $\pm 5.35$  mmHg

Tolerance in cuff measurement:  $\pm 3$  mmHg

Tolerance error in pulse rate measurement:  $\pm 5\%$

Technical specifications subject to change without prior notification.

## Technical testing procedures

---

### A) Function testing

Function testing of the unit can only be performed on human subjects or using a suitable simulator.

### B) Leak testing of the pressure circuit and error testing of the display

#### Notes:

- a) For testing to be carried out, the tube connecting the cuff must be pulled off from the connector plug.

The long connection nozzle on the connector plug, previously attached to the connecting tube, is now turned to connect it to the cuff socket on the unit. Make sure that the connector plug is pushed in as far as it can go i.e., that the seal on the connection nozzle is inside the cuff socket.

The remaining connection nozzle on the connector plug, previously attached to the unit, is connected either to the Normal (when testing the pressure reading for error) or to the cuff and a pump bulb with valve (when leak testing the pressure circuit).

- b) If the pressure accidentally exceeds 320 mmHg when performing the pressure test, an acoustic

signal is generated and the measurement value in the display blanks out every other second. Even if the pressure continues to increase, the pressure indicated remains at 320 mmHg.

### Testing

1. Switch off the unit.
2. Set the pressure preselector switch to 150 mmHg.
3. Press the START/STOP button **momentarily**.
4. When „0“ is flashing set the pressure selector switch to 240 mmHg and then back again to 150 mmHg. The unit is in test mode when the display indicates SYSTOLE, DIASTOLE and PULSE 0 mmHg.

5. Perform the usual tests for

- tolerance testing of the pressure reading.
- leak testing of the pressure circuit (take account of the cuff's settling time of one minute).

6. Switch off the unit and restore the original connections.

### **C) Security**

For security purposes the two halves of the unit can be joined together with a security mark. Alternatively the hole for the casing connector screw on the rear side of the unit can be secured by sealing it.

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**BOSCH + SOHN GmbH u. Co. KG**